

PS DETECT Order Form

PATIENT NAME _____ DOB _____

MRN (IF AVAILABLE) _____

PHYSICIAN NAME _____

PHYSICIAN ADDRESS _____


PHYSICIAN SIGNATURE _____ DATE _____

PLEASE SELECT ONE:

- ☐ PS DETECT TEST
- ☐ POST VARICOCELE REPAIR
- ☒ SEMEN

SPECIMEN TYPE:

Patient Instructions



Scan QR code to order your PS Detect kit through the PS Fertility portal.*

PS Detect kit ships directly to your home in discreet packaging.

Collect your sample in the privacy of your own home and ship back to PS Fertility.

PS Detect results will be provided to you and your physician.

*If this is a follow up test, please contact support@psfertility.com to obtain a QR code for a discounted test.




LEARN MORE AT [PSFERTILITY.COM](https://psfertility.com)

CONTACT AT [SUPPORT@PSFERTILITY.COM](mailto:support@psfertility.com)

An Exciting New Way to Detect Fertility

PS Detect is a novel male fertility test that detects whether sperm have the ability to fertilize an egg.



LOW PS → LOW FERTILIZATION RATE
Even if a standard semen analysis is normal.

SEMEN ANALYSIS IS MORE COMPLETE WITH FERTILIZATION COMPETENCY
PS Detect results can help guide treatment options.